









S No	PRODUCT ID	PRODUCT NAME	Image	PRODUCT FULL DESCRIPTION
1	GeIBSS	Gel Grip Exerciser/Stress & Strain Reliever (Round Ball Type)-SOFT-SMALL		<p>Muscular stress, strain & pain in the joints are associated with work hazard and infrequent injuries that occur in our regular life. It is recommended to use regularly our ergonomic designed Gel Grip Exerciser. One can minimize the harmful effects of stress on muscles and joints. The product effectively releases the stress, strain and pain of the type, neuro- muscular, forearm and palm disorders and post operative surgical or fracture in the forearm including phalange.</p> <p>Usage:</p> <ol style="list-style-type: none"> 1. Place the soft Gel Grip Exerciser in full open palm and hold with a gentle grip. 2. Compress and release gently to the convenience and comfort of user, as many a times daily, regularly. 3. Gel Grip Exerciser also tends to relieve stress and strains of working executives by keeping them more relaxed after the strenuous days work. 4. This product is washed with a mild soap, rinsed under normal water and shade dry. 5. It do not absorb liquid and take bad odor. <p>Note:</p> <ol style="list-style-type: none"> a. Keep off sharp/pointed implements. It does not bear to heavy objects. Do not expose to heat or direct sunlight. Do not stretch. Can be damaged if you tug it
2	GeIBSL	Gel Grip Exerciser/Stress & Strain Reliever (Round Ball Type)-SOFT-LARGE		
3	GeIBHS	Gel Grip Exerciser/Stress & Strain Reliever (Round Ball Type)-HARD-SMALL		
4	GeIBHL	Gel Grip Exerciser/Stress & Strain Reliever (Round Ball Type)-HARD-LARGE		
5	GeIGSS	Gel Grip Exerciser/Stress & Strain Reliever (kidney shape)-SOFT-Standard		






6	GeIINS	Full Gel Insole pair (Size: Small, Medium, Large)		<p>RECOMMENDED FOR: Gives comfort to the foot sole area in advancing age. The increase in body weight effects the foot by uneven distribution of body weight leading to severe foot pain. In the foot deformity it is readily available for its correction. Doctors have essential advice for post operative foot sores for cushion comfort and even distribution of the body weight.</p> <p>BENEFITS: The pressure exerted on the foot sole area is absorbed by the Ezy Walk Gel Insole and evenly distributes it, further do not allow to generate any pressure points leading to pain, permitting the person pain free walking with cushioning comfort.</p>
7	GeIHC	Gel Heel Cup pair (Size: Small, Large)		<p>RECOMMENDED FOR: Heel bears the maximum weight of the body while standing. Plantar Fascia produces uncanny pain which is a very common orthopaedic problem. The inevitable part of the treatment is to provide heel cushioning in the inner area of the shoe. Normal footwear do not have this facility. The development of the calcaneum spur and the overgrowth which is un-operatable, presses hard on the innersole heel area resulting unbearable pain and restricts the mobility.</p> <p>BENEFITS: Ezy-Walk Heel Cup are designed to provide four directional cushioning effect which regains its shape and size (do not deform) even after prolonged usage, thus imparts natural and effective cushioning interface between calcaneum and inner shoe base.</p>



8	GelAS	Gel Arc Support pair (Size: Small, Medium, Large)		<p>RECOMMENDED FOR: The flat foot is too common problem in which the person cannot take a straight posture. The absence of the mid-half cave-in do not allow to put a total support on the foot area, therefore the inability to stand erect without pain and strain is a common phenomenon located in the person's with this defect.</p> <p>BENEFITS: The Ezy-Walk arc support is a readymade solution. Usually footwear's do not have this facility. Preparing a "Whole Footwear" with inbuilt arc support is time consuming and not so cost effective. Ezy-Walk arc supports which is Gel based, non-toxic & non-allergic can be affixed easily to any existing footwear in use. This treatment is highly recommended by the physiotherapist and doctors.</p>
9	GelHP	Gel Heel Pad (Size: Standard)		<p>RECOMMENDED FOR: Heel bears the maximum weight of the body while standing. Plantar Fascia produces uncanny pain which is a very common orthopaedic problem. The inevitable part of the treatment is to provide heel cushioning in the inner area of the shoe. Normal footwear do not have this facility. The development of the calcaneum spur and the overgrowth which is un-operatable, presses hard on the innersole heel area resulting unbearable pain and restricts the mobility.</p> <p>BENEFITS: Ezy-Walk Heel Cup are designed to provide four directional cushioning effect which regains its shape and size (do not deform) even after prolonged usage, thus imparts natural and effective cushioning interface between calcaneum and inner shoe base.</p>



10	GelToe	Gel Toe Spreader / Bunion Cushion		<p>Placed between the first two toes, these Soft Gel Toe Spreaders straighten and align the big toe with the big toe joint to aid in the treatment and prevention of painful bunions. These TOE SPREADERS can also be used to space out and straighten other misaligned toes of the feet-and can help relieve the pain of overlapping toes by separating them.</p>
11	GelSeat12	Gel Seat Cushion 12" x 12" Eco (0.5" Plain Gel on top of 1" foam base with full lycra cover)		<p>Fatal injuries under severe shocks, coma etc., where mobility restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury. Soft comfortable, durable, seat/sofa/chair cushions. Having high density foam base with Gel layering covered with stretchable cloth .</p> <p>Recommended for computer seat cushion and extensive long time sitting for personnels in offices. Reduces the pressure and strain on the back bone and gives a soothing effects to the hips and thighs enabling tension free working.</p>
12	GelSeat14	Gel Seat Cushion 14" x 14" Eco (0.5" Plain Gel on top of 1" foam base with full lycra cover)		<p>Fatal injuries under severe shocks, coma etc., where mobility restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury. Soft comfortable, durable, seat/sofa/chair cushions. Having high density foam base with Gel layering covered with stretchable cloth .</p> <p>Recommended for computer seat cushion and extensive long time sitting for personnels in offices. Reduces the pressure and strain on the back bone and gives a soothing effects to the hips and thighs enabling tension free working.</p>






13	GelSeat16	Gel Seat Cushion 16" x 16" Eco (0.5" Plain Gel on top of 1" foam base with full lycra cover)		<p>Fatal injuries under severe shocks, coma etc., where mobility restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury. Soft comfortable, durable, seat/sofa/chair cushions. Having high density foam base with Gel layering covered with stretchable cloth .</p> <p>Recommended for computer seat cushion and extensive long time sitting for personnels in offices. Reduces the pressure and strain on the back bone and gives a soothing effects to the hips and thighs enabling tension free working.</p>
14	GelSeatR16	Gel Seat Cushion 16" x 16" Eco (1" Ridged/Pipe Type Gel on top of 1" foam base with full lycra cover)		<p>Fatal injuries under severe shocks, coma etc., where mobility restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury. Soft comfortable, durable, seat/sofa/chair cushions. Having high density foam base with Gel layering covered with stretchable cloth .</p> <p>Recommended for computer seat cushion and extensive long time sitting for personnels in offices. Reduces the pressure and strain on the back bone and gives a soothing effects to the hips and thighs enabling tension free working</p>
15	GelSeatBE16	Gel Seat Cushion 16" x 16" ELASTOMERIC Ball type seat cushion with elastomeric covering		<p>Fatal injuries under severe shocks, coma etc., where mobility restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury. Soft comfortable, durable, seat/sofa/chair cushions. Having high density foam base with Gel layering covered with stretchable cloth .</p> <p>Recommended for computer seat cushion and extensive long time sitting for personnels in offices. Reduces the pressure and strain on the back bone and gives a soothing effects to the hips and thighs enabling tension free working</p>






16	GelCoccxHoles	Gel Coccyx Seat Cushion – Holes all over (17" x 15" x 1.5" x 3.5")		<p>The Coccyx Cushion is ideal for prolonged use as a chair cushion in the office, as a car seat cushion or at home. Not only does it relieve pressure on the Coccyx (tail bone), but it improves sitting posture by raising the hips in relation to the knees. Improved sitting posture can relieve and prevent low back pain when sitting.</p> <p>The Coccyx Cushion is especially suitable for people who have a bruised or painful tail bone region. It is also appropriate for the treatment and prevention of lower back pain caused by disc problems and poor sitting posture.</p>
17	GelCoccxR	Gel Coccyx Seat Cushion – Pipe Type (17" x 15" x 1.5" x 3.5")		
18	GelCoccxB	Gel Coccyx Seat Cushion – Ball type (78 Gel Balls on Top) (17" x 15" x 1.5" x 3.5")		






19	GelToilet	Gel Commode Seat Pad		<p>Ideal padded commode seat cushion that can be used on nearly all bathroom commode safety products. The Gel cushioning provides pressure relief while toileting. Ideal for patients who use toilets, bed pans or commodes for extended periods.</p>
20	GelIOTComp	<p>Operation Table GEL Mattress (Set of 7 pcs.) a. 1 pc. (19.5" x 18"), b. 1 pc. (19.5" x 18"), c. 1 pc. (19.5" x 17"), d. 1 pc. (19.5" x 4"), e. 1 pc. (19.5" x 4"), f. 1 pc. (19.5" x 4"), g. 1 pc. (14" x 9")</p>		<p>All incidences of operations have painful and stress causing effect on the patient. The procedure involves lot of time during labor and major surgical cases causing strain to the patient. The comfortable Gel cushioning provides intense relief and cooperation by the patient to carry out the functional prospects of operation.</p>
21	GelBaby	Gel Baby Pillow		<p>Every baby can benefit from Gel Support that's affordable and Long-Lasting. Gel pillows are an effective way to minimize pressure and provide uniform support for the soft cranium of the pre-term infant. The wonderfully soft surface and supple support of Gel Pillows keeps babies more comfortable and minimizes the re-positioning needed to alleviate pressure points. Can be used for micro preemies, preemies, and full-term babies ... even as they grow.</p>






22	GelSheet12	Soft Gel Sheet 12" x 12" x 0.12" (with fabric lining)		<p>Fatal injuries under severe shocks, coma etc., where bed restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury.</p> <p>Soft, comfortable and durable soft Gel Sheetings. Having inbuilt fabric lining for better adhesion.</p> <p>Recommended for layering in various orthopaedic and rehabilitation aids where anti vibration, shear and cushioning are required</p>
23	GelSheet1224	Soft Gel Sheet 12" x 24" x 0.12" (with fabric lining)		<p>Fatal injuries under severe shocks, coma etc., where bed restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury.</p> <p>Soft, comfortable and durable soft Gel Sheetings. Having inbuilt fabric lining for better adhesion.</p> <p>Recommended for layering in various orthopaedic and rehabilitation aids where anti vibration, shear and cushioning are required</p>
24	GelSheet2424	Soft Gel Sheet 24" x 24" x 0.12" (with fabric lining)		<p>Fatal injuries under severe shocks, coma etc., where bed restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury.</p> <p>Soft, comfortable and durable soft Gel Sheetings. Having inbuilt fabric lining for better adhesion.</p> <p>Recommended for layering in various orthopaedic and rehabilitation aids where anti vibration, shear and cushioning are required</p>






25	GelSheet2436	Soft Gel Sheet 24" x 36" x 0.12" (with fabric lining)		<p>Fatal injuries under severe shocks, coma etc., where bed restrictions is compulsion for long period, Soft Gel cushioning prevent developing pressure points which might add to injury.</p> <p>Soft, comfortable and durable soft Gel Sheetings. Having inbuilt fabric lining for better adhesion.</p> <p>Recommended for layering in various orthopaedic and rehabilitation aids where anti vibration, shear and cushioning are required</p>
26	GP-Neuro	Gel Base for Metallic Neuro Head Frame (pair)		<p>Conforms comfortably to the patient's head in metallic neuro head frame. Standard design, can be attached to any existing frame. Distributes weight & relieves pressure effectively. Special plain Gel-without elastomeric coating</p>
27	GP-OHR	Open Head Ring (Horse Shoe type) (200 x 75 x 45) mm		<p>Soft, comfortable and durable, the Gel Open Head Ring supports the patient's face, while allowing easy access for intubation during surgery.</p>






28	GP-CHRS	Gel Positioners - Closed Head Ring, Adult (200 x 75 x 45) mm		<p>Conforms comfortably to the patient's head. Distributes weight & relieves pressure effectively. Suitable for use in many surgical positions.</p>
29	GP-CHRL	Gel Positioners - Closed Head Ring, Adolescent (140 x 60 x 35) mm		<p>Conforms comfortably to the patient's head. Distributes weight & relieves pressure effectively. Suitable for use in many surgical positions.</p>
30	GP-Sac	Gel Positioners - Sacral Protector (440 x 445 x 30) mm		<p>The Sacral area is a high risk pressure area if left unprotected. The sacral protector has been designed to provide support in the supine and lithotomy positions by relieving pressure and focusing on peripheral support.</p>






31	GP-HPL	Gel Positioners - Hip Pad (520 x 520 x 10) mm - Large		<p>Used to protect the body during procedures in the Supine, Prone, Trendelburg, Reverse Tendelburg and Lithotomy positions.</p>
32	GP-HPS	Gel Positioners - Hip Pad (400 x 400 x 10) mm - Small		<p>Used to protect the body during procedures in the Supine, Prone, Trendelburg, Reverse Tendelburg and Lithotomy positions.</p>
33	GP-UBNPL	Gel Positioners - Ulnar/Brachial Nerve Protector (450 x 200 x 10) mm - Large		<p>Provides excellent support and stability to ulnar/brachial nerves during surgery.</p>





34	GP-UBNPS	<p>Gel Positioners – Ulnar/Brachial Nerve Protector (350 x 150 x 10) mm – Small</p>		<p>Provides excellent support and stability to ulnar/brachial nerves during surgery.</p>
35	GP-HSPL	<p>Gel Positioners – Heel Support Pads (200 x 100 x 65) mm (Pair) Large</p>		<p>Used to protect the body during procedures in the Supine, Prone, Trendelburg, Reverse Tendelburg and Lithotomy positions.</p>
36	GP-HSPS	<p>Gel Positioners – Heel Support Pads (150 x 75 x 50) mm (Pair) Small</p>		<p>Used to protect the body during procedures in the Supine, Prone, Trendelburg, Reverse Tendelburg and Lithotomy positions.</p>



37	GP-FBCRL	Gel Positioners - Flat Bottom Chest Roll (510 x 150 x 125) mm – Large		Provides excellent support and stability during surgery in the prone position.
38	GP-FBCRS	Gel Positioners - Flat Bottom Chest Roll (510 x 100 x 75) mm – Medium		Provides excellent support and stability during surgery in the prone position.
39	GP-PPM	Gel Positioners - Patient Positioner Medium (295 x 110 x 60) mm		Provides excellent support and stability during surgery in the prone and other positions.



40	GP-PPL	Gel Positioners - Patient Positioner Large (380 x 140 x 80) mm		Provides excellent support and stability during surgery in the prone and other positions.
41	Pros-Breast	Gel External Breast Prosthesis		These inexpensive realistic breast forms for mastectomy are an alternative to the more costly silicone breast forms / breast prostheses or other mastectomy products. They have most of the features of the silicone mastectomy products with the EXCEPTION OF THE COST. These prosthetic breast forms are produced in standard and plus sizes as per the requirement and can be easily fitted in the existing bra-cup.

Terms:

- Prices are subject to change without any prior notice
- Before placing order kindly confirm the dispatch lead time

Manufactured by:

GEL O' KARE

39, Adarsh Vihar-II, Lane Adjacent to Nobel Public School, 1 Km. from Dubagga Taxi Stand, Mahipatmau, Hardoi Road, Lucknow-227 017 Uttar Pradesh, INDIA
 Phone: +91 522 3258325, 098397 58992, 094150 89959
 Web Site: www.gel-o-kare.com E-mail: gelokare@gmail.com, vermaashwani@hotmail.com